

### **Western Blot Detection of Protein**

1. Obtain from your TA a strip of western blot nitrocellulose membrane. The TA placed the strip into a 15 ml tube and added 5 ml of blocking solution (buffer 1) to the blot.
2. **Start here!** carefully pour off buffer 1, discarding the solution. (Note: all solutions are discarded after they are used, but save the blot!)
3. Add 2 ml (1:5,000) of mouse anti-V5 epitope primary antibody (**1°**) onto the blot and incubate with shaking for 30 minutes. Make sure the blot is laying flat against the tube and is covered with solution. Tape the tube on a shaking platform.
4. Pour off the antibody solution, add 5 ml of buffer 2, shake the blot for 3 minutes. Pour off the wash solution and repeat the washing step one more time. (These steps will wash off unbound antibodies).
5. Pour off the last wash solution, add 2 ml of anti-mouse alkaline phosphate conjugated secondary antibody (**2°**) (1:2500). Incubate at room temperature with the secondary antibody for 15 minutes.
6. Pour off the previous solution and incubate the nitrocellulose filter strip with 5 ml of buffer 2 at room temperature for 3 minutes. Repeat this wash one more time.
7. Pour off the previous solution and incubate the nitrocellulose filter with 5 ml of Buffer 3 at room temperature for 1 minute (this step will equilibrate the filter to the proper pH for hybrid detection).
8. Pour off the solution, add 2 ml of the color detection reagent. Let the color develop under a sheet of foil for 10 minutes to overnight, do not shake the blot. (The length of incubation time depends on the abundance of the protein on the blot).

**Buffer 1: (Blocking buffer and antibody binding buffer)**

<u>Reagents</u>	<u>Final Concentration</u>
Maleic acid	0.1M
NaCl	0.15M
BSA	1%
Triton X-100	0.3%
pH 7.5	

**Buffer 2: (TBS Tween, washing buffer)**

<u>Reagents</u>	<u>Final Concentration</u>
Tris -HCl	0.1M
NaCl	0.15M
Tween 20	0.1%
pH 7.5	

**Buffer 3: Alkaline phosphatase detection buffer**

<u>Reagents</u>	<u>Final Concentration</u>
Tris -HCl	0.1M
NaCl	0.1M
MgCl <sub>2</sub>	0.05M
pH 9.5	

**Color detection reagent**

2 ml of buffer 3 with BCIP and NBT (Promega).